OPTIKA

M I C R O S C O P E S

I T A L Y

Ver. 2.0.0

ryf ag

Ryf AG
Bettlachstrasse 2
2540 Grenchen
tel 032 654 21 00
fax 032 654 21 09

www.ryfag.ch

B-600 MET

 **OPERATION MANUAL**

 **GUIDA UTENTE**

 **MANUAL DE INSTRUCCIONES**

 **MANUEL D'INSTRUCTIONS**

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

1.0 DESCRIPTION	page 3
2.0 INTRODUCTION	page 5
3.0 UNPACKING AND ASSEMBLY	page 5
4.0 USING THE MICROSCOPE	page 7
5.0 MAINTENANCE	page 10
6.0 ELECTRICS	page 11
7.0 RECOVERY AND RECYCLING	page 12

Metallographic features

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily classroom and laboratory use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

3.0 UNPACKING AND ASSEMBLY

- 3.1** The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.
- 3.2** The objectives are located in individual protective vials. Remove the objectives from the vials and insert them into the microscope nosepiece in the order from the lowest magnification to the highest, in a clockwise direction from the rear.
- 3.3** In a separate box you find the metallographic attachment. Take it out and place it onto the top of the microscope's stand. Then lock it in position by turning the screw with the provided allen wrench:

Locking screw for the attachment

3.0 UNPACKING AND ASSEMBLY

- 3.4 Place the observation head into the top of the attachment and tighten the lock-screw.
- 3.5 Insert the eyepieces into the eye tubes.
- 3.6 Insert polarizer and analyzer filters into their slots in the metallographic attachment.
- 3.7 Connect the power cable of the epillumination to the appropriate connector on the rear:

Connector
for LED
epillumination

4.1 Adjust the observation head

Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw.

4.2 Place the specimen on the stage

Put the specimen on the mechanical stage (transparent glass plate). Ensure that the specimen is centred over the stage opening by adjusting the coaxial knobs of the stage.

Brightfield vision

4.3 Lamp settings

The microscope is fitted with a 50 W illuminator. Insert the plug of the cable into the power socket and turn on the switch on the back of the main body. Turn the brightness adjustment knob (on the right of the body) until you hear a “click”, meaning that the lamp is on, then keep turning the knob to a brightness level suitable for observation.

Lamp housing also has centring knobs for fine adjustment of illumination.

4.4 Adjust interpupillary distance

Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. After finishing the interpupillary adjustment, turn the dioptic adjustment ring on the left eye tube to the “0” position at the graduation line on the eyetube.

4.5 Focus and diopter adjustment

Loosen the focus-stop knob (near the left focus knob), rotate the coarse focusing knob to bring the slide into focus with an objective with a low magnification, and then lock the focus-stop knob. Adjust the fine focusing knob to get the image sharp and clear while observing with your right eye, then turn the left dioptic adjustment ring to a sharp and clear image also with the other eye. Turn the tension-adjust-knob (near the right focus knob) to get a suitable tension for the focus system. When the image is in focus, turn the nosepiece to choose the objective you need.

4.6 Condenser adjustment

Turn the condenser adjustment knob to move the condenser up or down for a clear and uniform sample illumination. The two condenser-centring screws can be used to centre the condenser.

4.7 Set the numerical aperture

Adjust the aperture of the iris diaphragm to set the numerical aperture of the illuminator, thus controlling image contrast and resolution.

4.8 Field diaphragm

Adjust the field diaphragm to align the system and to control stray light.

4.9 Additional filters

The chromatic correction filter (blue) can be inserted into the holder (just above the ring of the field diaphragm) for definition of specimen parts.

4.10 Video capturing (optional)

B-600MET be connected to cameras via a photo/video adaptor, for photo and video capturing. Before taking a picture or filming video, if necessary, obscure with a dark cloth both the camera viewfinder and the eyepieces and pull out the light path selector lever. Please refer to the adaptor and camera manuals for further details.

Epillumination vision (metallography)

4.11 Lamp settings

The metallographic attachment is fitted with an high power white LED. Insert the plug of the cable into the power socket and turn on the switch on the back of the main body. Turn the brightness adjustment knob (on the left of the body) until you hear a “click”, meaning that the LED is on, then keep turning the knob to a brightness level suitable for observation. Please note: if no light comes out from the objective, check that both diaphragms on the metallographic attachment are open (levers UP).

4.12 Adjust interpupillary distance


Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. After finishing the interpupillary adjustment, turn the dioptic adjustment ring on the left eye tube to the “0” position at the graduation line on the eyetube.

4.13 Focus and diopter adjustment

Loosen the focus-stop knob (near the left focus knob), rotate the coarse focusing knob to bring the slide into focus with an objective with a low magnification, and then lock the focus-stop knob. Adjust the fine focusing knob to get the image sharp and clear while observing with your right eye, then turn the left dioptic adjustment ring to a sharp and clear image also with the other eye. Turn the tension-adjust-knob (near the right focus knob) to get a suitable tension for the focus system. When the image is in focus, turn the nosepiece to choose the objective you need.

4.14 Field diaphragm

Adjust the field diaphragm to control stray light. If you want to close the diaphragm, pull down the lever:


At the bottom of the lever, on both sides of the microscopes, there is a couple of screws that can centre the field diaphragm.

4.15 Set the numerical aperture

Adjust the aperture of the iris diaphragm to set the numerical aperture of the epilluminator, thus controlling image contrast and resolution. If you want to close the diaphragm, pull down the lever:

At the bottom of the lever, on both sides of the microscopes, there is a couple of screws that can centre the aperture diaphragm.

4.16 Polarization analysis

Slightly pull outwards the polarizer filter (marked "S-PO") to insert it into the optical path. Turn the rotating disc of the analyzer filter: if the two filters are correctly inserted, you should find a position of the analyzer where no light is transmitted ("crossed polarizers" position).

4.17 Darkfield vision

After you focus the sample with the proper objective, slowly push the darkfield lever toward the inside. Stop pushing as soon as you see a shadow covering your sample. In this manner you should obtain an oblique illumination that greatly enhances the reliefs and superficial defects of your sample.

5.1 Always think about

- The following environment is required: Indoor temperature: 0-40°C, Maximum relative humidity: 85 % (non condensing).
- Keep the microscope away from dust and shocks while in use.
- Turn off the light immediately after use.
- Use a soft lens tissue to clean the optics after use.
- Only if needed, use a cloth moistened with water and a mild detergent, rinsing with water and drying immediately with a lint-free cloth.
- After use, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

5.2 Do not!

- Wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Use solvents, neither on the microscope, nor on the optics.
- Disassemble objective or eyepieces to attempt to clean them.
- Mishandle or impose unnecessary force on the microscope.
- Clean the unit with volatile solvents or abrasive cleaners.
- Attempt to service the microscope yourself.

5.3 How to change the brightfield lamp

- Before changing the lamp, you must pull the plug out of the electrical socket. To avoid being burnt, wait until the lamp has cooled down.
- Remove the lamp housing by pulling it out of the back of the microscope.
- Pull out the old lamp.
- Hold the new lamp without touching it with your bare hands (e.g. using its packaging, or a cloth) and insert it into the holder. Replace/close the lamp holder.


- 5.4** If you need to send the microscope to Optika for maintenance, please use the original packaging.

6.0 ELECTRICS

Mains Power:	90-240 V, 50/60 Hz
Brightfield Lamp:	12V/50W halogen lamp
Epillumination:	High power white LED

The appliance reports the symbol:

This symbol means that the appliance can be a precious source of raw materials. Therefore, it must not be disposed of as waste, but separately collected for the recycling and the recovery of the materials it contained in it. Such materials, if improperly dispersed into the environment, can be harmful to the environment and to human health.

The producer of the equipment, Optika Microscopes, recovers, re-uses and recycles the raw materials contained in the equipment. Such recovery, however, needs your help.

When, at the end of its operating life, you shall decide to dispose of the apparatus, do not try to open it, nor to use parts of it in ways other than reported in this Manual, but bring it back to the Distributor you bought it from, or to other Optika Microscopes distributors. The Distributor shall collect the apparatus free of charge.

The recovery of the raw materials shall then take place in accordance with the European Directive 2002 / 96 / EC and all other relevant Directives. Never disassemble, nor dispose of as waste, apparatus reporting the “crossed bin” symbol indicated above.

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