## RLQ1101

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Vielen Dank

# RLQ1101

This high quality product has been designed and developed to require a minimum of maintenance. Please follow the guidelines below and read this manual in its entirety to ensure safe operation of the unit.

## Warning

- Please always pull mains plug before cleaning the device.
- Storing and using the cold light source only in dry rooms.
- The provided power cord has a grounded plug. Please always use power cable with correctly earthed socket.

## Safety Information

- To avoid eye damage, don't look directly into the light source when turned on.
- To avoid burns, please let the lamp cool down before replacing a spare part.

## Symbols

The following symbols are used in this manual.

**CAUTION!** This symbol indicates a potential risk and alerts you to proceed with caution.

## Intended Use

The RLY1101 Cold Light Source is perfectly designed for most of the technical applications where you need professional and exact lighting solutions, e.g. in quality management, laboratories, industrial endoscopy, vision systems and many more.

The light source impresses by a modern and state-of-the-art design, which stands out by its modularity and its flexibility. It can be positioned in a horizontal and vertical way.

The RLQ1101 is a high quality light source, which provides the user a perfect light quality for an attractive price.

## Package Contents

Description	Quantity
Cold Light RLQ1101	1
Power Cord	1
Empty Filter Holder	1
User Manual	1

## Getting Started

The cold light source is a plug & play device very easy to use.


Please check before use whether ...

- ... the mains plug is connected.
- ... the main switch at the back side is switched on.  
(blue LED at the front must glow)
- ... the brightness is set to the desired value.
- ... the areas for the air intake (front) and out take (rear) are free.
- ... the location has a good ventilation in order to avoid a thermal overload of the system.

## Description of Buttons and Switches


### Front View (fig.1)

- a) **Electronic dimming**  
For continuous dimming of the light and for the complete switch off of the lamp.
- b) **Mechanical dimming with iris diaphragm**  
For continuous dimming of the light without change of colour temperature.
- c) **Light guide access**  
Connection of the light guide with a standard diameter of 15mm.
- d) **Air intake slots**  
Air cooling must be ensured at any time.
- e) **Screw caps**
- f) **Cover of housing**


### Back View (fig.2)

- g) **Plug for power cord**
- h) **Fine fuse**  
The fine fuse can be removed by pulling out the insert with a small tool
- i) **Main switch**  
For switching on and off of all functions of the light source.
- j) **High quality low-noise rear fan**  
The air out take must be ensured at any time

### Interior (fig.3)

- k) **Cover plate**
- l) **Connecting clamp**
- m) **Lamp socket**
- n) **Halogen projection lamp**
- o) **Lamp holder**

# RLQ1101

## Troubleshooting

The information in the table below can help to resolve the possible operating problems quickly.

Problem	Causes	Solution
No function at the device. No blue LED:	No network connection.	Please check, whether the power cord is connected.
	Fine fuse defective.	Please check the fine fuse (see picture). Replace the fuse in case of a defect. Please only use fuses with the right electrical values.
No light, but fan in operation. Blue LED is turned on:	Soiled air intakes could cause a too high temperature of the device. Automatic switch off by automatic thermo protection.	Please ensure a secure air intake. Clean the air intake and outtake areas.
	A too high environment temperature.	The environment temperature should not exceed 40° C max.
	Projection lamp is defective.	Please change lamp according to user manual.

## Storage, Maintenance, Repair

The device is maintenance-free. Some points, however, you should take into consideration.

1. Please unplug the device before opening and cleaning
2. For the cleaning of the housing please use a smooth, semi-wet cloth, if necessary with a smooth household cleaning fluid
3. No moisture should enter into the interior of the housing.
4. Before replacing the spare lamp please clean the air intake and outtake areas with a small brush.
5. To avoid damages during the storing and transport please retain the packaging
6. Repairs are only allowed by trained and experienced personnel. In case of inappropriate use the warranty expires.

## Technical Service

### Web Resources:

Visit our Ryf website at [www.ryfag.ch](http://www.ryfag.ch) for:

- Contact details of technical customer service.
- Ryf product overview and price list as well as infos about accessories and spare parts.
- Additional product information and special offers.

### Contact:

For information or technical assistance contact:  
Ryf AG, Bettlachstrasse 2, CH-2540 Grenchen.

## Replacement of Lamps

1. Please switch off the halogen projection lamp with the electronic dimming switch (a).
2. Let the cold light source cool down for several minutes by using the fan.
3. Switch off the cold light source with the main switch (i).
4. Disconnect the power plug (g).
5. Unscrew the four screw caps at the surface of the device (e).
6. Remove the cover of the housing (f).
7. Open the cover plate (k).  
**Warning!** The lamp could be hot. **Danger of burns!**
8. Ensure that the projection lamp has been cooled down.
9. Remove the projection lamp from the lamp holder (o).
10. Remove the lamp socket (m).
11. Please only use special halogen projection lamps 100W/12V (Art.No.100-000365) for replacement. **Attention!** The use of other lamps could cause damage in your device.
12. Clean the air intake and outtake areas (d and j) with a smooth tissue when replacing the lamp.
13. Please touch the new lamp only at the outside of the lamp.
14. Put the lamp socket (m) onto the new projection lamp.
15. Put the projection lamp in the lamp holder (o).
16. Close the inside cover plate (k).
17. Close the housing with the cover (f).
18. Fix the screw caps on the housing (e).
19. Plug in the power cord (g).
20. Now, your cold light source is ready to use.

**Further repairs are only allowed by trained and experienced personnel. In case of improper handling the warranty expires.**

Technical Specifications	
Weight:	2,3kg
Dimensions:	260 x 175 x 135mm ( L x W x H )
Operating voltage:	100-240VAC / 50/60Hz
Halogen lamp:	100W / 12V
Power consumption:	111W
Average lamp life span:	1.500h
Light guide access:	standard 15mm
Light control:	electronic and mechanic continuous adjustment
Colour temperature:	3.100K
Material:	aluminium
Colour:	nature and titan anodized
Cooling system:	active cooling by a low-noise-fan

Accessories and Spare Parts	
	<ul style="list-style-type: none"> <li>Daylight filter - Art.-No.: 100-003962</li> <li>Green filter - Art.-No.: 100-003960</li> <li>Yellow filter - Art.-No.: 100-003958</li> <li>Red filter - Art.-No.: 100-003959</li> <li>Blue filter - Art.-No.: 100-003961</li> </ul>
	<p>Wide range of standard light guides and individual customized solutions.</p>
	<ul style="list-style-type: none"> <li>Fuse 250V / 1,25A - Art.-No.: 100-001484</li> <li>Lamp socket GZ 6,35 - Art.-No.: 100-002749</li> </ul>
	<p>Special halogen projection lamp 100W / 12V - Art.-No.: 100-000365</p>

## Warranty

Ryf AG warrants that this product will be free from defects in material and workmanship for a period of one (1) year from date of purchase. If a defect is present, Ryf AG will, at its option, repair, replace, or refund the purchase price of this product at no charge to you, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear.

For your protection, items being returned must be insured against possible damage or loss. This warranty shall be limited to the replacement of defective products.

It is expressly agreed that this warranty will be in lieu of all warranties of fitness and in lieu of the warranty of merchantability.

## Equipment Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste. Instead, it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you

## Notizen/Notes

### Ryf AG

Bettlachstrasse 2, CH-2540 Grenchen, Switzerland  
Tel.: +41-(0)32-654 2100, Fax: +41-(0)32-654 2109  
E-Mail: ryfag@ryfag.ch, Website: www.ryfag.ch