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MODELLE

RLQ 1100 LED RLQ 2600 LED

STAND AUGUST 2012

I INHALTSVERZEICHNIS

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1 DEVICE OVERVIEW

FRONT VIEW RLQ 1100 LED (SIMILAR TO RLQ 2600 LED)

1	Main switch (ON/OFF)	CHAP. 3.3
2	Incremental encoder	CHAP. 3.4
3	LCD display	CHAP. 3.5
4	Light guide holder with collet chuck and safety shutdown	CHAP. 3.1
5	Filter slider	CHAP. 3.6
6	Carrying handle	
7	Ventilation openings (device rear)	
8	Ventilation openings (device bottom)	
9	Nameplate (device bottom)	
10	Insert nuts for optional assembly (device bottom)	

**DEVICE REAR
ZLED CLS 6000**

**DEVICE REAR
ZLED CLS 9000**

- 7 Ventilation openings
- 11 USB interface
- 12 Footswitch interface
- 13 RS232 interface (only with ZLED CLS 9000!)
- 14 CAN interface (only with ZLED CLS 9000!)
- 15 Connection socket for power cable

2 IMPORTANT NOTES

SYMBOLS:

SYMBOL	DESCRIPTION
	WARNING OF A HAZARD, OBSERVE DOCUMENTATION!

USE TO THE INTENDED PURPOSE:

The cold light source RLQ 1100 LED / RLQ 2600 LED is intended for indoor use in industry and laboratories. The cold light source must not be used in potentially explosive rooms or environments.

When used to the intended purpose, the cold light source serves to illuminate objects of all kinds by directing high-intensity visible light through various light guides to the object.

The cold light source RLQ 1100 LED / RLQ 2600 LED meets the requirements of the following standards and directives:

- DIN EN 61010-1 (IEC 61010-1)
- DIN EN 62471
- DIN EN 61326-1
- Directive 2004/108/EC
- Directive 2006/95/EC

Compliance of the device with requirements in the above-mentioned standards and directives is evidenced by means of the pertinent technical documentation.

SAFETY NOTES:

CAUTION! The cold light source emits visible high-intensity light. Measurement of photo-biological safety according to DIN EN 62471 resulted in the classification of the device into Risk Group 2.

Please study and observe these Operating Instructions carefully, as they contain important safety information relating to the safety of the user and the device. Safety of the device and the user is not ensured if the Operating Instructions are not adhered to.

CAUTION! Potentially hazardous optical radiation. Avoid looking into the lamp or the light guide aperture emission for longer periods during operation. May be harmful to the eyes.

Light-absorbing materials convert impinging light into heat. Hence, heat-sensitive and/or inflammable light-absorbing materials may suffer thermal damage and/or burns. Please observe the below notes in order to avoid such damage and hazards:

- Do not cover open light guide holders or light guide outlets at any time! Risk of fire!
- Avoid at all cost covering open light guide holders or light guide outlets by hands or other parts of the body. Risk of burns!
- When illuminating heat-sensitive or inflammable light-absorbing objects (e.g. in microscopy), make sure that there is a distance between the relevant material and the employed light guides, and adjust the intensity of the light source so as to prevent damage and/or hazards to the object.
- All light guide outlets that are not in use during operation must be kept at a safe distance – i.e. at least 10 cm – from heat-sensitive or inflammable light-absorbing materials while the light source is switched on.
- Biological tissue may be exposed to unnecessarily high stress and may be subject to damage during excessive illumination with visible light, hence the intensity of the light source and the duration of illumination of such tissues should be limited to the absolutely necessary levels and periods.

The cold-light source may only be operated at the voltage indicated on the nameplate (9).

Ensure a safe protective conductor connection.

The cold-light source has been developed for indoor use in dry rooms.

The cold-light source must not be used in potentially explosives environments.

Safe disconnection from power supply can only be ensured by disconnecting the mains plug.

All ventilation openings (7, 8) must be uncovered at all times to ensure sufficient cooling of the device. In case of insufficient cooling an integrated sensor will continuously control light intensity down until shutdown to protect the illuminant from overheating and damage. In this case, the LCD (3) will issue an error message (see Chapter 3.5. for LCD messages).

The device may only be opened or disassembled by authorised service personnel. Technical changes on the device are prohibited. Repair may only be effected by manufacturer or authorised service personnel.

These Operating Instructions must be swiftly available and accessible to any device operator.

Manufacturer does not assume liability for any damage that results from failure to observe these Operating Instructions.

NOTES FOR WASTE MANAGEMENT:

Please note that you should not dispose of the cold light sources RLO 1100 LED and RLO 2600 LED with domestic waste. If you would like to dispose of your device please send it back to Ryf AG (address See Reverse site). Ryf AG will dispose of your device professionally and harmlessly to the environment.

3 OPERATION

3.1 LIGHT GUIDE SOCKET

Rotate the collet chuck of the light guide holders (4) in counter clockwise direction to open the collet chuck. Push the light guide into the light guide holder to its limit-stop and rotate the collet chuck (4) in clockwise direction to close it and to ensure that the light guide is firmly positioned. Light guides with guide pin: Make sure to insert the guide pin into one of the guiding slots.

CAUTION! The light guide holder (4) is equipped with a safety shutdown to protect the user from dazzling. If the light guide is inserted into the light guide holder incorrectly or not all the way to the limit-stop, the light of the light source will not start when switched on, and the message 'Light Guide!' is shown in the display.

3.2 MAINS CONNECTION

A power cable is not included in delivery. A country specific power cable is available from authorized dealers of the light source or directly from ZETT OPTICS GmbH.

Push the barrel connector of the power cable into the connection socket (15) at the device rear provided for this purpose. Connect the plug connector of the cold device cable to power mains.

3.3 START-UP

Switch on the cold-light source by actuating the main switch (1). The LCD display (3) of the cold light source is active in switched-on condition. Upon each switch-on, the light source starts with the intensity setting of the prior use.

CAUTION! Pull the mains plug for safe disconnection of the device from the power supply!

3.4 INCREMENTAL ENCODER

Rotate the incremental encoder (2) in order to adjust intensity (in clockwise direction: up; in counter-clockwise direction: down).

By repeatedly briefly pressing the incremental encoder (2) the intensity can be adjusted in pre-set stages of 0, 20, 40, 60, 80, and 100% (factory pre-settings).

A menu can be accessed by pressing the incremental encoder (2) for at least 3 seconds. The menu navigation is shown in the LCD display (3). In the menu, individual intensities can be set and stored in 6 memory locations. In addition, fine setting of the incremental encoder can be determined by the user, to be chosen from the settings 1 (100 steps) up to 5 (20 steps).

The factory settings of the device can be re-instated as well. The menu functionality is explained in detail in Chapter 3.5 "LCD display".

3.5 LCD DISPLAY

Menu navigation for RLQ 1100 LED, functions are also applicable to RLQ 2600 LED.

3.5.1 MENU NAVIGATION

```
RLQ 1100
Intens : 50%
```

Intensity has been set to a value (e.g. 50%). Rotate the incremental encoder (2) to change the intensity. Press the incremental encoder repeatedly to set the intensity values 0, 20, 40, 60, 80, and 100% quickly.

```
RLQ 1100
Intens : 33%
```

Intensity has been set to 33%. Rotate the incremental encoder (2) to vary intensity. Press the incremental encoder (2) to access the menu.

```
*** MENU ***
SET M1 >
```

By rotating the incremental encoder (2) you get to the settings 1 to 6 („SET M1-6“) Pre-settings are 0%, 20%, 40%, 60%, 80%, and 100%.

```
*** MENU ***
SET M2 >
```

Briefly pressing the incremental encoder (2) again accesses a sub-menu, where new intensity values can be set in the 6 memory locations Set M1 to Set M6 by rotating the incremental encoder (2).

The new value (e.g. 5%) can be saved by briefly pressing the encoder.

```
*** MENU ***
<SET M2 5
```

```
*** MENU ***
DimSteps >
```

By rotating the incremental encoder (2) you get to the menu item „DimSteps“. Pre-setting of „DimSteps“ is 2. That means that the intensity can be varied in 2% steps by rotating the incremental encoder (2).

```
*** MENU ***
<DimSteps:1%
```

Briefly press the incremental encoder (2) to access a sub-menu where the fineness of intensity stages can be adjusted. Here you can choose from 1% up to 5%. Confirm the selection by briefly pressing the encoder again

```
*** MENU ***
Reset all >
```

Rotate the incremental encoder (2) in the menu to access the menu item 'Reset all'.

```
Reset all?
< NO
```

Briefly press the incremental encoder (2) to access a sub-menu where the factory settings of the device can be re-instated. The display shows the question 'Reset all?'. Rotate the incremental encoder to decide between 'NO' and 'YES' and briefly press the incremental encoder (2) to confirm the selection.

```
Reset all?
< YES
```

```
*** MENU ***
FW-Info >
```

By rotating the incremental encoder (2) you get to the menu item „FW-Info“.

```
Firmware
Ver.: X.XXX
```

Briefly press the incremental encoder (2) to get informationen about softwareversion („Firmware Ver. xxx) and serialnumber of the device.

```
Serialnumber
```

```
*** MENU ***
< Exit Menu
```

Rotate the incremental encoder (2) in the menu to access the menu item 'Exit Menu'. Briefly press the incremental encoder (2) to exit the menu and revert to normal operation mode.

3.5.2 OTHER MESSAGES

RLQ 1100
* Standby *

With the footswitch (accessory) connected the light source can be switched on and off quickly. After switch-off of the light source by the footswitch, 'Standby' will be shown in the display.

3.5.3 ERROR MESSAGES

RLQ 1100
Light Guide!

If no light guide was inserted or it was not inserted up to the limit-stop, the light of the light source does not start upon switch-on and the display (3) shows the message 'Light Guide!'. In such case, please, insert a suitable light guide.

RLQ 1100
Temp-Sensor!

The integrated temperature sensor is not functioning properly, i.e. there is no protection from potential overheating of the device. Switch off the device and contact our customer service department.

RLQ 1100
LED open!

The LED ist defective or power supply is interrupted. Switch off the device and contact our customer service department.

RLQ 1100
LED Overtemp

The illuminant has heated up too much. The illuminant was switched off automatically. Switch off the device and contact our customer service department.

RLQ 1100
PCB Overtemp

The control PCB has heated up too much. The illuminant was automatically switched off. Switch off the device and contact our customer service department.

RLQ 1100
Fan error!

The fan has become blocked or is defective. Sufficient cooling is not ensured. Switch off the device and contact our customer service department.

3.6 FILTER SLIDER

The filter slider (5) is fitted with 3 positions, the first two of which can be pulled out and equipped with the appropriate colour and conversion filters designed for this purpose. The third position is an empty position that cannot be equipped with a filter. The filter slider (5) slightly locks into place when the respective position has been reached.

CAUTION! Only use original filters for RLQ 1100 LED / RLQ 2600 LED and ensure that the filter slider lock into the desired position in order to avoid damage. The light source emits high-intensity light. Use of unsuitable filter may cause thermal damage of the filters or the device. Thermal damage may be caused on the filter slider if the filter is not locked into one of the three positions.

Original filters can be bought from authorized dealers of the light source or directly from Ryf AG.

3.7 INTERFACES

RLQ 1100 LED

The RLQ 1100 LED is fitted with an USB interface (11) and an interface to connect a foot-actuated switch (12). The interfaces are located at the device rear and are marked with the lettering 'USB' and 'Foot Switch'.

The light source can be connected via the USB interface to a PC with a USB interface. You can find software for demonstration on the added CD. With this software you can adjust the intensity of the cold light source. Users who are aware of programming have the opportunity to reprogram the software for their needs.

A footswitch can be connected to the footswitch interface in order to switch the light of the light source on and off. A suitable footswitch is available from authorized dealers of the light source or directly from ZETT OPTICS GmbH.

RLQ 2600 LED

The RLQ 2600 LED is fitted with an USB interface (11), a CAN interface (14), a RS232 interface (15), and an interface to connect a footswitch (12). The interfaces are located at the device rear and are marked with the lettering 'USB', 'CAN', 'RS232', and 'Foot Switch'.

The light source can be connected via one of the interfaces to a PC with the same interface. You can find software for demonstration on the added CD. With this software you can adjust the intensity of the cold light source.

Users who are aware of programming have the opportunity to reprogram the software for their needs.

A footswitch can be connected to the footswitch interface in order to switch the light of the light source on and off. A suitable footswitch is available from authorized dealers of the light source or directly from Ryf AG.

3.8 ACCESSORIES

A comprehensive accessories programme for RLQ 1100 LED / RLQ 2600 LED is available from authorized dealers of the light source or directly from Ryf AG.

Faultless function, safety and light yield of the light source are only ensured when original Light guides by Ryf AG are used.

CAUTION! Use of other light guides may trigger the safety cut-out which prevents switching on the illuminant of the cold-light source. In such case the message 'Light Guide!' is shown in the display die.

4 MAINTENANCE

The ZLED CLS 6000/ ZLED CLS 9000 is a maintenance-free device.

The cold light source is not designed for disinfection (e.g. in medical applications). Please clean the cold light source only externally with a soft dry cloth. Slightly moisten the cloth with a commercial plastic cleaning cloth to remove tenacious soiling.

CAUTION! Do not spray any fluids into the ventilation openings (7, 8). Do not use a wet cloth to clean the cold-light source.

5 TROUBLESHOOTING

TYPE FAULT	POSSIBLE CAUSE	REMEDIAL ACTION
Illuminant fails to light up	No light guide inserted. Unsuitable light guide inserted. Light guide not inserted into light guide holder deep enough to the limit-stop	Insert a suitable light guide into the light guide holder
	Device not switched on	Switch on the device
	Plug not in socket	Put the plug into the socket outlet
	No voltage	Check voltage
Luminous intensity deteriorates	Electronic parts overheated	Ensure sufficient cooling, restart the device after a longer cooling-down period

In case of any malfunctions of your cold light source that cannot be rectified by means of the above measures, please contact your specialised dealer or Ryf AG (address see Reverse site).

CAUTION! Do not open the device yourself. Repairs may only be carried out by authorised and skilled personnel.

6 TECHNICAL DATA

ZLED MODEL CLS 6000

GENERAL INFORMATION

Type designation	RLQ 1100 LED
Dimensions (W x H x D) [mm]	192 x 111 x 236
Weight [g]	3095
Cooling	Axial-flow fan
Temperature: in normal operating mode	5-40° C
Temperature during transport/storage	-40-70° C
Air humidity, pressure in normal operating mode	10-92%, no condensation
Air humidity, pressure during operation	30-75%
Atmospheric pressure in normal operating mode	760-1060 hPa
Atmospheric pressure during transport/storage	700-1060 hPa
Use to intended purpose	Use in interior rooms

ELECTRICAL DATA

Operating voltage, frequency	100–240V ± 10%, 50–60 Hz
Overvoltage category	II
Power consumption max.	70 W
Degree of protection	IP 20
Degree of soiling	2

LIGHTING DATA

Maximum usable active diameter of light guide [mm]	9
Total luminous flux at light guide outlet (ZETT OPTICS light guide flexible, active Ø 9 mm, L 1000 mm), stage 100, [lm]	600

EMC EMISSION CLASS

B

ZLED MODEL CLS 9000**GENERAL INFORMATION**

Type designation	RLQ 2600 LED
Dimensions (W x H x D) [mm]	192 x 111 x 236
Weight [g]	3095
Cooling	Axial-flow fan
Temperature: in normal operating mode	5-40° C
Temperature during transport/storage	-40-70° C
Air humidity, pressure in normal operating mode	10-92%, no condensation
Air humidity, pressure during operation	30-75%
Atmospheric pressure in normal operating mode	760-1060 hPa
Atmospheric pressure during transport/storage	700-1060 hPa
Use to intended purpose	Use in interior rooms

ELECTRICAL DATA

Operating voltage, frequency	100–240V ± 10%, 50–60 Hz
Overvoltage category	II
Power consumption max.	80 W
Degree of protection	IP 20
Degree of soiling	2

LIGHTING DATA

Maximum usable active diameter of light guide [mm]	9
Total luminous flux at light guide outlet (ZETT OPTICS light guide flexible, active Ø 9 mm, L 1000 mm), stage 100, [lm]	870

EMC EMISSION CLASS**B**

7 EXPOSURE HAZARD VALUES (EHV)

RLQ 1100 LED	With light guide	Without light guide
Photochemical hazard: Free group	41.3	60.7
Photochemical hazard: Risk group 1	5.9	4.3
Photochemical hazard: Risk group 2	0.01	0.01
Thermal hazard: Free group/ Risk group 1	2.5	1.8
Thermal hazard: Risk group 2	0.8	0.6

RLQ 2600 LED	With light guide	Without light guide
Photochemical hazard: Free group	52.5	88.6
Photochemical hazard: Risk group 1	8.8	5.9
Photochemical hazard: Risk group 2	0.02	0.01
Thermal hazard: Free group/ Risk group 1	3.7	2.5
Thermal hazard: Risk group 2	0.94	0.95

8 HAZARD DISTANCES (HD)

RLQ 1100 LED	With light guide	Without light guide
Photochemical hazard: Free group	1.3m	1.6m
Photochemical hazard: Risk group 1	0.5m	0.4m
Thermal hazard: Free group/ Risk group 1	0.3m	0.3m

RLQ 2600 LED	With light guide	Without light guide
Photochemical hazard: Free group	1.5m	1.9m
Photochemical hazard: Risk group 1	0.6m	0.5m
Thermal hazard: Free group/ Risk group 1	0.4m	0.3m

4 FOR YOUR NOTES

RYF AG

Bettlachstraße 2
2540 Grenchen / Switzerland

phone: +41 32 654 21 00
fax: +41 32 654 21 09

info@ryfag.ch
www.ryfag.ch | www.fotoryf.ch | eshop.fotoryf.ch