**RLQLED01**  
**LED-Kaltlichtquelle**  
**LED Cold Light Source**

**BEDIENUNGSANLEITUNG**  
**USER MANUAL**

**Artikelnummer / Article number**  
**100-004987**

**Version: 1.0**

## RLQLED01

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Vielen Dank

## RLQLED01

This high quality product has been designed and developed to require a minimum of maintenance. Please follow the guidelines below and read this manual in its entirety to ensure safe operation of the unit.

**Caution!** This symbol indicates a potential risk and alerts you to proceed with caution

### Warning

- Please disconnect plug from the main socket before cleaning the device.
- For indoor use only

### Safety Information

- To avoid eye damage, don't look directly into the light source during operation

## Intended Use

The LED cold light source RLQLED01 represents for almost all usual applications a cost efficient and energy saving alternative to conventional halogen light sources. Contrary to halogen cold light sources it transfers a much larger part of the energy to visible light and has a life span of about 50.000 hours. At the light guide the RLQLED01 couples in roughly as much as light as a 150W halogen light source, however only consumes about 20W.

The RLQLED01 is perfectly designed for most of the technical applications where you need professional and exact lighting solutions, e.g. in quality management, laboratories, industrial endoscopy, machine vision systems and many more.

The light source impresses by a modern and state-of-the-art design which stands out by its modularity and its flexibility. It can be positioned in a horizontal and vertical way and can be easily adapted to individual OEM versions.

Package Contents	Quantity
Cold Light Source RLQLED01	1
Power Supply Plug, 24V with „industry plug“	1
User Manual	1

## Getting Started

The Roma LED-xv is a plug & play device very easy to use.

Please check before use whether ...

- the power cord is connected
- the device is switched on (the blue LED has to glow).
- the brightness is set to the desired value


## RLQLED01

### Description of Buttons and Switches

#### Front View (fig.1)

- a) **Elektronic, digital dimming**  
For the dimming of the light (standby until 100%) and for the activation of the memory function
- b) **Light Guide Access**  
Connection of the light guide with a standard diameter of 15mm
- c) **Operation Control**


(fig.1)

#### Rear View (fig.2)

- d) **Connection to the power adaptor**

(fig.2)

### Digital Potentiometer with Memory Function

The brightness can be adapted with a digital potentiometer. The highest level is indicated by a short flickering of the LED.

Memory function:

1. On the switch off of the device with the potentiometer the last setting is automatically saved for the next use.
2. In case of an interruption of the power supply (e.g. withdrawing of the power adaptor or switch off of the main switch) the last setting can be saved individually by the user:
  - At the preferred level please push the potentiometer for 2,5s.
  - The light goes off for this period and appears again after 2,5s
  - The desired setting is saved.
  - The power supply can be interrupted.

## RLQLED01

### Troubleshooting

Review the information in the table below to troubleshooting operating

Problem	Causes	Solution
Power supply:	No network connection.	Please check whether the power adaptor is connected.
No function at the device. No blue LED.	Electronic defective.	Please contact your distributor.
No light, blue LED glows.	Dimming function on stand by.	Please adapt light to the desired value.
	LED and/or electronic defective.	Please contact your distributor.
LED gets dark or yellowly	Thermal overload of the system.	Please ensure an ambient temperature of 0°-40° C.
	LED life span exceeded.	Please contact your distributor.

### Repair and Maintenance

The device is basically free of maintenance. Some points, however, you should take into consideration.

1. Please unplug the device before cleaning
2. For the cleaning of the housing please use a smooth, semi-wet cloth, if necessary with a smooth household cleaning fluid
3. No moisture should get into the interior of the housing.

**Further repairs are only allowed by trained and experienced personnel. In case of improper handling the warranty expires.**

### Technical Service

#### Web Resources:

Visit the StarLight Website on [www.starlight-sl.de](http://www.starlight-sl.de) for:

- Contact information for the technical service
- StarLight product overview and price list as well as additional product information and special offers.

#### Contact:

For information or technical assistance contact your local distributor or directly StarLight Sensorik Lichtleitsysteme GmbH & Co. KG.

## RLQLED01

Technical Specifications	
Weight:	3,45kg
Dimensions:	250 x 174 x 81mm (L x W x H)
Operating voltage:	24V DC, protection class III (safety extra-low voltageSELV)
LED colour:	White, similar to day light
Colour temperature:	Typically 5.600K
Ø LED life span:	About 50.000 h (according to manufacturer)
LED:	15W power LED, white
Light guide access:	Standard 15mm
Colour reproduction index:	80
Active diameter of the light guide:	max. 8mm
Light control:	Electronic, digital dimming
Material:	Aluminium
Operating temperature:	0°C - 40°C
Relative humidity::	max. 95%, non condensated
Cooling system:	Passive cooling without fan
Colour of housing:	Natur und titan anodised

Accessories	
	Power adaptor plug for UK Art. No.: 100-005047
	Wide range of standard light guides and individual, customized solutions.

Spare Parts	
	Power adaptor plug 24V with industry plug. Operating voltage: 100-240V, 50-60Hz Output voltage: 24V DC Protection rating II (protection insolated) Art. No.: 100-004900

## RLQLED01

### Warranty

StarLight Sensorik Lichtleitsysteme GmbH & Co. KG warrants that this product will be free from defects in material and workmanship for a period of two years from date of purchase. If a defect is present, StarLight will, at its option, repair, replace, or refund the purchase price of this product at no charge to you, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear.

For your protection, items being returned must be insured against possible damage or loss. This warranty shall be limited to the replacement of defective products. It is expressly agreed that this warranty will be in lieu of all warranties of fitness and in lieu of the warranty of merchantability.

### Equipment Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste. Instead it's your responsibility to correctly dispose of your equipment at lifecycle -end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you

## Notizen/Notes

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---

- 

---