



Inverted Material microscope Set IMM-1 (30W)

OPTIKA

M I C R O S C O P E S
I T A L Y



INDUSTRY SERIES

Günstiges Inverses Materialmikroskop-Set mit 30 Watt Beleuchtung
Microscope inversé métallographique avec éclairage à 30 Watt




Inverted Material microscope Set IMM-1 (30W)

IMM-1

NEW

With infinitive optical system, providing excellent optical functions.
 With the features of stable stand structure, advanced stage design and comfortable operation, it is widely used in institute and laboratory to observe and identify the structure of various metal and alloy.

Viewing head	18-75 mm TYrinoocular head inclined 30°, interpupillary 18-75mm.
Eyepiece	EW10X/20 high point, extra wide field eyepiece EW10X/20
Objective	4X/0.1 10X/0.25 20X/0.40 40X/0.60
Nosepiece	Quintuple nosepiece
Stage	Plain stage 160x250mm
Focusing	Coaxial coarse and fine adjustment, vertical objective movement coarse stroke: 37.7 mm for rotation, fine stroke: 0.2mm for rotation.
Kohler illumination	Halogen lamp 6V30W
Filter	Blue, yellow, green and ground glass



Inverted metallurgical microscope