



**Model:**  
CL-15.1

**Typology:**  
ILLUMINATOR

**Description:**  
*LED ring light illuminator.*

**Cromaticity coordinates (x,y):** 0.320, 0.320 (pure white), 6300K

**Optical power:** Luminous flux: 220 lumen  
Illuminance (at 10 cm distance): 30000 lux

**Power supply:** Buttons for brightness adjustment.  
Different illumination settings: omnidirectional (ring completely on), oblique (only a quarter or half ring switched on).  
External power supply provided.  
Input: 110/240 Vac 50/60 Hz; output: 12 Vdc 500 mA

**Accessories:** Fixing ring for stereo heads.

ryf ag



Ryf AG  
Bettlachstrasse 2  
2540 Grenchen  
tel 032 654 21 00  
fax 032 654 21 09

[www.ryfag.ch](http://www.ryfag.ch)