

Universal LED transmitted light stand for stereo microscopes

Ryeco RSZ-1012

ryeco
microscopes
+ swiss quality



Ryeco tripod RSZ-1012 with ring lighting Ryf NKL-12 LED

Universal LED transmitted light stand for stereo microscopes



Ryeco universal transmitted light stand with very good price-performance ratio:

Inexpensive transmitted light stand for various stereo microscopes such as Ryeco RSZ-45(T), Nikon SMZ 445 / SMZ 745(T), Leica S9, Zeiss 305 / 508 and Olympus SZ51 / SZ61(T).

The two built-in lights can be operated ergonomically using the two buttons in the base plate:

- LED transmitted light with glass plate, adjustable
- Side LED adjustable, angle of the LED can be tilted
- Main switch can be operated from the right
- Optionally with ring lighting Ryf NKL12-LED / RL4b

The stand with integrated Z-axis drive is made of aluminium and is painted white (origin Asia).

Dimensions: W=245mm x D=285mm x H=340mm

Ryeco Swiss Warranty & Service Package & Final Assembly

SAP Order number: RSZ-1012

www.ryfag.ch