

Universal Ryeco stand set for smaller stereo microscopes

Ryeco RSZ-STL1H

ryeco
microscopes
+ swiss quality



Ryeco simple universal stand RSZ-STLH with the ring light Ryf NKL-12 LED
(The rack and pinion drive head holder must be ordered separately with no. 235520)

Universal Ryeco stand set for smaller stereo microscopes

Ryeco compact universal stand with head holder / drive with very good price-performance ratio:

Affordable entry-level universal stand for various stereo microscopes such as Ryeco RSZ-45(T), Nikon SMZ 445 / SMZ 745(T), Leica S9, Zeiss 305 / 508 and Olympus SZ51 / SZ61(T).

This simple industrial-grade Universal tripod, in a set with focussing drive 235120, has no built-in lighting, but can be supplemented with the most common LED ring light illuminators such as the Ryf RL4b, NKL-12 LED or the most cost-effective variant, the CL-16.1.

Depending on the model, an RGL-4 or an adapter is also required.

The base is made of steel and painted. The two column poles are also made of steel and are chrome-plated. The entire tripod is sold as a set with a head holder with rack and pinion drive.

The entire drive is made of aluminium and is painted white / grey (both from Asia).



Dimensions: Base plate = 230x230 / Pillar Ø 32mmx380mm / Horizontal arm = 430mm

Ryeco Swiss Warranty & Service Package & Final Assembly

SAP Order number: RSZ-STL1H

SAP Order number: 235520